



2020 Product Catalog

ALL ENFit® ALL the TIME!





Vesco

Latin verb meaning "to feed"

Feeding Innovation

Our promise to deliver a continuous stream of break-through
enternal products to better serve our customers





High-Quality Tube-Feeding Devices

- ISO 9001 & 13485 Certified Quality Management System
- FDA Registered & 510(k) Premarket Clearance
- Made without natural rubber latex, DEHP or other phthalates plasticizers

Unbeatable Service

- Friendly, skilled and experienced customer care specialists
- Full portfolio of ENFit products direct or through distributors
- Responsive when and where you need us

Competitive Pricing

- Lean, agile organization driven to keep costs down
- High-quality products without the premium price
- Volume discounts for even more savings


ENFit Tip Syringes

ENFit Tip Syringes for Feed, Flush & Irrigation (5, 10, 20, 35 & 60mL)


- ISO 80369-3 compliant ENFit connector designed to reduce the risk of tubing misconnections
- Enteral specific syringes for use in initial tube placement, flushing/hydration, bolus feeding and irrigation
- Printed graduated length markings in mL
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers
- One cap included in each blister pack

ENFit Low Dose Tip Syringes (1 & 3 mL)


- Designed to minimize process changes in pharmacy and bedside administration
- Demonstrated dose accuracy administering smaller volumes compliant with GEDSA*
- Internal lumen consistent with traditional male syringe design reduces deadspace and priming volume
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers
- Sterile blister packs with one cap included in each pack



ENFit Standard Tip



ENFit Low Dose Tip (LDT)



ENFit Syringe with Transition Connector

| THUMB RING IRRIGATION ENFIT SYRINGES | Item # | Units/Carton |
|--------------------------------------|---------|--------------|
| 60mL | VED-661 | 180 |
| 60mL + 500mL Container | VED-665 | 180 |

ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission.

* GEDSA August, 2016 Guidance Document

Syringe Item Numbers


| ENFIT TIP SYRINGES | | | CLEAN (Non-Sterile) | | |
|--------------------|-------------------|------------|---------------------|---------------------------|--------------|
| Size (mL) | STERILE | | with Cap | with Transition Connector | Units/Carton |
| | Item # | Units/Case | | | |
| 1 (LDT) | VED-601EO w/ Caps | 800 | N/A | N/A | N/A |
| 3 (LDT) | VED-603EO w/ Caps | 800 | N/A | N/A | N/A |
| 5 | VED-605EO | 800 | VED-605 | VED-605TC | 500 |
| 10 | VED-610EO | 800 | VED-610 | VED-610TC | 500 |
| 20 | N/A | N/A | VED-620 | VED-620TC | 500 |
| 35 | VED-635EO | 400 | VED-635 | VED-635TC | 500 |
| 60 | VED-660EO | 400 | VED-660 | VED-660TC | 180 |



ENFit Pharmacy Items

ENFit Medication Bottle Fill Caps (Multiple Sizes)

- Pharmacy friendly purple cap designed to aid in single handed placement
- Various sizes to meet most common medication bottles
- Fits inside bottle opening to allow continued use of child proof caps
- Ideal for aseptic technique
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers



ENFit Transition Connectors

- Transition Connector consistent with the design and taper of a catheter tip syringe
- Available prepackaged within the inner case or within each syringe package
- Also available as a stand alone item by the case

Transition Connectors

ENFit Syringe Tip Cap (Fits All Sizes)

- Pharmacy friendly purple cap designed to aid in single handed placement
- Ideal for aseptic technique and secure for transport
- Made without natural rubber latex, DEHP, Bisphenol A or other phthalates plasticizers



ENFit Medicine Straw (50 mm)

- Connects to an ENFit Tip Syringe to easily facilitate the withdrawal of medication or feed

Important Safety Information

- Use ENFit tip syringes only with ENFit feeding tube ports. Feeding tubes with different open female ports are not compatible with Vesco ENFit tip syringes
- ENFit tip syringes have the potential to misconnect with small bore connectors of other healthcare applications
- Single-patient use and Rx only



ENFit Medicine Straw

ENFit Medication
Bottle Adapters