

Adapt Medical Adhesive Spray

General Information

- Medical grade silicone adhesive
- Increases the adhesion of skin barriers and/or tape
- Helps to provide a secure seal

Using Medical Adhesive Spray

1. Prepare your pouching system as instructed.
2. Remove the protective release liner from the adhesive (skin barrier and/or tape) on your pouching system.
3. Shake **Adapt** medical adhesive spray can well before using.
4. Spray directly on to the adhesive side of the pouching system. Allow to dry for 5 minutes.
5. Apply your pouching system as instructed.

Frequently Asked Questions

Q: Is it necessary to use the Adapt medical adhesive spray on the entire skin barrier and tape?

A: You can spray it on the entire adhesive surface of the pouching system. However, if there are particular areas that are not adhering well, you can use it only in those areas.

Q: What is the best way to spray the adhesive?

A: Hold can so that the arrow faces the product to be sprayed. Spray in short bursts until the surface of the product is coated evenly.

Q: Can a second coat be applied?

A: Yes, if necessary. Remember to wait 5 minutes for the second coat to dry before applying your pouching system.

Q: Can I spray it directly on my skin?

A: This product should not be sprayed directly on skin and it is not necessary in order for it to be effective.

Q: Will the adhesive be difficult to remove from the skin?

A: The medical adhesive should come off easily when the pouching system is removed. If needed, the Adapt medical adhesive remover spray or Adapt universal remover wipe can help with ease of removal.



Adapt Medical Adhesive Spray
Stock No 7730



Spray directly on to the adhesive, and allow to dry for 5 minutes. Then apply pouch as instructed.

Adapt Medical Adhesive Remover Spray/ Adapt Universal Remover Wipe

General Information

- Helps remove adhesive skin barriers and/or tape
- Removes adhesive residue from the skin
- Silicone-based formulation
- No sting from alcohol



Adapt Medical Adhesive Remover Spray/Adapt Universal Remover Wipe
Stock No 7737/7760

Using Adhesive Remover Spray

1. Shake well before using.
2. Peel back one edge of the skin barrier or tape border.
3. Hold the can so that the arrow faces the base of the skin barrier or tape, and spray until the adhesive releases. The 360° spray technology minimizes waste, and will allow you to spray from any angle.
4. Ensure your skin is dry and apply your new ostomy pouching system as usual.



Spray
Spray until the skin barrier releases.

Using Universal Remover Wipes

1. Tear open the package and take out the wipe.
2. Press or rub the wipe against your skin as you peel off your skin barrier and/or tape.
3. Press or rub the wipe over any areas of the skin where tape residue remains.
4. Ensure your skin is dry and apply your new ostomy pouching system as usual.



Wipes
Rub the wipe over skin as you peel off skin barrier.

Frequently Asked Questions

Q: Is an adhesive remover always needed?

A: No, this product is optional. Use it to remove adhesive residue from your skin, or to make it easier to remove your skin barrier or tape.

Q: Will this be uncomfortable to use?

A: This silicone-based product is formulated to be a no-sting alternative. When used correctly, the adhesive remover will help to make the removal of your pouching system easier and help reduce the potential for skin stripping.



Wipes
Clean off remaining tape residue.

Q: Do I need to wash the adhesive remover off of my skin?

A: The adhesive remover evaporates quickly and dries without residue, so the next barrier can attach securely. If you wash your skin after use be sure to rinse and dry thoroughly before the next pouch application.

For product questions, sampling needs, or detailed clinical questions concerning our products in the US, call **1.888.740.8999**. In Canada call **1.800.263.7400**.

Routine follow-up with your healthcare professional is recommended.

Prior to using any ostomy products/accessories be sure to read all product inserts and labels for any indications, contraindications, hazards, warnings, cautions and directions for use.