

THE OBVIOUS CHOICE

People with swallowing problems prefer SimplyThick® over other thickeners. The reasons are simple: no graininess, no aftertaste, and better tasting beverages.* In addition, beverages look like they should.

SimplyThick is extremely easy and convenient to use. Any time, anywhere, with any beverage, SimplyThick will work for you. No messy powder. No bulky cans. No measuring spoons.

Try the samples you have received and see why patients prefer SimplyThick. Then, see the enclosed details on our special 50% rebate offer and purchase SimplyThick for yourself!



DISPENSER BOTTLE USAGE CHART

Desired Consistency	Per 4 oz./120 ml	Per 8 oz./240 ml	Per 32 oz./960 ml
Nectar	1 Full Stroke	2 Full Strokes	8 Full Strokes
Honey	2 Full Strokes	4 Full Strokes	16 Full Strokes
Pudding	4 Full Strokes	8 Full Strokes	32 Full Strokes

PREPARING SIMPLYTHICK

SimplyThick's pre-measured packets allow for fast and easy preparation of beverages anywhere. Each packet produces either nectar or honey consistency. **Two honey packets produce pudding consistency.**

The SimplyThick® brand thickener dispenser bottle is designed for simple, quick use in situations that require maximum flexibility in thickening. Simply add the appropriate amount of gel and mix.

Individual preparation

15g & 30g Individual (smaller) packets thicken 4 oz./120 ml

step 1

Measure 4 oz./120 ml of beverage for each packet of SimplyThick, or dispense SimplyThick into beverage according to the chart on the left.

mixing

Mix briskly with fork-like beating an egg for 20-30 seconds.

step 2

shaking

Shake vigorously for 5-10 seconds, then enjoy SimplyThick!

mixing

Mix briskly with whisk-like beating an egg for 20-30 seconds.

step 3

Ensure there is at least 1-2" (2.5-5cm) of head space in container and put lid on container.

Ensure there is at least 2-3" (5-7.5 cm) of head space.

Shaking hot or carbonated beverages is NOT recommended.

Mixing SimplyThick with a spoon is NOT recommended.

Bulk preparation

120g & 240g Bulk (larger) packets thicken 32 oz./960 ml

step 1

Measure 32 oz./960 ml of beverage for each packet of SimplyThick, or dispense SimplyThick into beverage according to the chart on the left.

mixing

Ensure there is at least 2-3" (5-7.5 cm) of head space.

step 2

Ensure there is at least 2-3" (5-7.5 cm) of head space.

step 3

Mix briskly with whisk-like beating an egg for 20-30 seconds.

Blenders and hand mixers are also great tools for bulk mixing SimplyThick.

SIMPLYTHICK® 50% REBATE

Up to a maximum rebate of \$20.00

To receive your 50% rebate** on SimplyThick:

1. Talk with your physician or healthcare professional to see whether SimplyThick is right for you.
2. Purchase any size of SimplyThick.
3. Complete this coupon and mail it with proof-of-purchase (UPC code from box and a copy of your receipt) to:

SimplyThick Rebate Offer
200 South Hanley Road
Suite 1102
St. Louis, MO 63105

Offer limited to one rebate per address.



[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!