

# SPECIMEN COLLECTION



DYND30330

## 4 oz. Basic Specimen Container

» 4 oz. (120 ml)

Item No.	Description	Pkg
DYND30330	Sterile Path, Basic Specimen Container, 4 oz.	100 ea/cs
DYND30331	O.R. Sterile, Basic Specimen Container, 4 oz.	100 ea/cs
DYND30335	Non-Sterile, Basic Specimen Container, 4 oz.	500 ea/cs



DYND30350

## 4.5 oz. Basic Specimen Container

» 4.5 oz. (130 ml)

Item No.	Description	Pkg
DYND30350	Sterile Path, Basic Specimen Container, 4.5 oz.	100 ea/cs
DYND30351	O.R. Sterile, Basic Specimen Container, 4.5 oz.	100 ea/cs



DYND30362

## 3 oz. Pneumatic Tube Safe Specimen Container

» Pneumatic tube safe with 95kPa specifications  
» 3 oz. (90 ml)

Item No.	Description	Pkg
DYND30362	Sterile Path, Pneumatic Tube Specimen Container, 3 oz.	400 ea/cs
DYND30367	O.R. Sterile, Pneumatic Tube Specimen Container, 3 oz.	100 ea/cs



DYND30364

## 4 oz. Pneumatic Tube Safe Specimen Container

» Pneumatic tube safe with 95kPa specifications  
» 4 oz. (120 ml)

Item No.	Description	Pkg
DYND30364	Sterile Path, Pneumatic Tube Specimen Container, 4 oz.	300 ea/cs
DYND30369	O.R. Sterile, Pneumatic Tube Specimen Container, 4 oz.	100 ea/cs



DYND30382

## 3 oz. Click & Close Specimen Container

» Secure closure is accomplished with a natural motion, and an audible “click” ensures a tight closure  
» Pneumatic tube safe with 95kPa specifications

Item No.	Description	Pkg
DYND30382	Sterile Path, Click & Close Specimen Container, 3 oz.	400 ea/cs
DYND30387	O.R. Sterile, Click & Close Specimen Container, 3 oz.	100 ea/cs



DYND30385

## 4 oz. Click & Close Specimen Container

» Secure closure is accomplished with a natural motion, and an audible “click” ensures a tight closure  
» Pneumatic tube safe with 95kPa specifications

Item No.	Description	Pkg
DYND30385	Sterile Path, Click & Close Specimen Container, 4 oz.	300 ea/cs
DYND30389	O.R. Sterile, Click & Close Specimen Container, 4 oz.	100 ea/cs